# WS-820

**INSTRUCTOIN MANUAL  
DIGITAL BLOOD PRESSURE MONITOR WS-820**

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
HUN

ROU

This manual contains information useful for safe and efficient operation of digital blood pressure monitor model WS-820 (hereinafter "Device"). Do not use this device in any manner other than that specified in this manual. Please look thoroughly through this manual and particularly "Recommendations for Correct Measurement".

## PARTS AND COMPONENTS

- 1 LCD DISPLAY
- 2 START BUTTON «» (START, ON/OFF)
- 3 BUTTON M1
- 4 BUTTON M2
- 5 BATTERY COMPARTMENT COVER
- 6 CUFF
- 7 CASE
- 8 BATTERIES

## GENERAL INFORMATION

### INDICATIONS FOR USE

WS-820 device is intended for measuring systolic and diastolic arterial pressure and for rating pulse in patients aged from 15 and elder with a cuff wrapped around a wrist. The device is recommended for use in patients with unstable (nonpermanent) blood pressure or hypertension at home as a supplement to medical surveillance. The cuff is suitable for a wrist with a circumference of about 12,5 - 21,5 cm. Blood pressure is measured in the range from 50 to 250 mmHg for systolic and 40 to 180 mmHg for diastolic, in the range from 40 to 160 heartbeats per minute.

The product is not designed for neonatal or infant use. Also, the effectiveness of the product is not established in person in pregnancy (including preeclamptic condition) as measurement accuracy has not been verified. Please consult with your doctor or physician to use this product for such subject group. The product is not intended for use during professional transport of a patient outside a healthcare facility.

## OPERATION PRINCIPLES

The device uses the oscillometric method of measurement. Connect a cuff to an electronic unit and wrap it around your wrist. When you press the «» button the device starts the automatic inflation, which is followed by the blood pressure measurement.

The sensor element of the device detects slight pressure oscillations in the cuff produced by the expansion and contraction of the brachial artery in response to each heartbeat. Rhythm and amplitude of the pressure waves are measured and displayed on the LC display as a numerical value in millimeters of mercury. The device has an irregular pulse rhythm indicator, as well as two memories with 30 cells in each calculating the average value.

## NISSEI New Technologies

**IQ SYSTEM** is an intelligent system developed by NISSEI engineers with due regard for the specific nature of each person's heartbeat, thus providing high accuracy of readings within minimal measurement time.

Irregular pulse rhythm indicator is a special icon on the display that informs on the irregular heartbeat, while the measurement result is correct. Please be noted that the irregular pulse can be caused by arrhythmia. To control this at an early stage we recommend you visit your doctor.

**M-Cuff** series is a unique cuff unit designed and patented by NISSEI. M-shaped cuff can read pulse wave simultaneously from two arteries.

**WS-820** is the most accurate instrument following up the results of the tests held by the German Hypertension Society\*. About 40 devices of different manufacturers were tested both on wrist and upper arm, and WS-820 got the highest score in terms of accuracy as of August 2006.

\* "GEK-Besser leben", 2/2006

## COMPLETE SET

The complete set WS-820 includes:

- Electronic unit with cuff – 1 pc.
- Batteries – 2 pc.
- Case – 1 pc.
- Instruction Manual – 1 pc.
- Warranty - 1 pc.
- Packaging – 1 pc.

## RECOMMENDATIONS ON CORRECT MEASUREMENTS

1 If treated with hemodialysis or anticoagulants, antiplatelets or steroids, refer to your doctor about the blood pressure measurement. Use of this instrument under such conditions could cause internal bleeding.

2 Malfunctions are possible when the device is used near working mobile phones, microwave ovens and other equipment generating electromagnetic radiation (e.g. The product should be used 3.3m away from the wireless communication device with a maximum output power of 2W).

3 Do not use this product along with implantable and wearable medical electrical equipment such as pacemaker, defibrillator, or electrocardiographic monitor. This product is not also intended to be used with HF surgical equipment as the use with such equipment may cause the malfunction of the product.

4 Do not use this product in an explosive environment such as near flammable anesthetics or inside oxygen chamber.

5 The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

6 Do not use any accessories other than those specified by the manufacturer. Otherwise, correct measurement readings cannot be obtained.

7 Do not apply the cuff over wounded wrist, wrist of which side is under an intravascular access or therapy or an arterio-venous shunt, or wrist on the side of a mastectomy. Otherwise injury may be resulted.

8 Make sure that inflation of the cuff is not causing prolonged impairment of blood circulation. Also, be cautious about temporary loss of the functions of any other medical equipment if any monitoring equipment is used on the same limb with the blood pressure measuring cuff.


9 Do not take out batteries when the monitor is turned on. Make sure that the monitor is turned off before removing batteries.

10 Do not inflate the cuff when it is not wrapped around your wrist.

11 This product is not designed for self-use by unspecified persons in public areas.

12 Do not apply the cuff on the limb which the intravenous drip infusion is implemented.

13 For correct measurement it is necessary to know that the BLOOD PRESSURE IS SUBJECT TO SHARP FLUCTUATIONS EVEN IN SHORT TIME INTERVALS. The blood pressure level depends on many factors. It is commonly lower in summer and higher in winter. Blood pressure varies along with atmospheric pressure and depends on the physical exertion, emotional excitability, stress and diet. Medical drugs, alcohol and smoking exert great influence as well. Occasionally, measurements in the clinic cause an increase in pressure values. Therefore, blood pressure measured at home is often different from that measured in the clinic. Since blood pressure increases at low temperatures, measurements should be made at room temperature (about 20°C). In case the product is stored in the environment with ambient temperature above 40°C or below 10°C, please leave it for at least 2 hours before taking a measurement.

During the day, the difference in the readings in healthy people may attain 30-50 mm Hg for systolic (upper) pressure and up to 10 mm Hg for diastolic (lower) pressure. Dependence of blood pressure on various factors is individual for each person. Therefore it is recommended to keep a special recording of blood pressure readings. **ONLY A DOCTOR MAY ANALYZE TRENDS IN CHANGING YOUR BLOOD PRESSURE BASED ON CORRESPONDING RECORDINGS.**


14 In case of cardiovascular diseases and a number of other diseases that require the blood pressure monitoring, measurements should be carried out in the hours specified by a doc-

tor. REMEMBER THAT THE DIAGNOSTICS AND ANY TREATMENT OF ARTERIAL HYPERTENSION SHOULD BE CARRIED OUT ONLY BY A DOCTOR BASED ON BLOOD PRESSURE READINGS OBTAINED BY A DOCTOR. MEDICAL DRUG ADMINISTRATION OR CHANGE OF DOSAGES SHOULD BE MADE ONLY BY PRESCRIPTION OF AN ATTENDING DOCTOR.

15 In case of disorders such as deep vascular sclerosis, weak pulse wave and break in rhythm of heart contractions, the correct blood pressure measurement can be complicated. IN THIS CASE, A DOCTOR SHALL PROVIDE RECOMMENDATIONS IN RELATION TO USE OF THIS DEVICE.

16 KEEP QUIET DURING THE MEASUREMENT TO OBTAIN THE CORRECT BLOOD PRESSURE READING WHEN USING THE ELECTRONIC DEVICE. The blood pressure measurement should be carried out in a quiet comfortable atmosphere at room temperature. Exclude meal an hour before the measurement, and exclude smoking, soft drinks, and alcohol 1.5-2 hours before the measurement.

17 Accuracy of the blood pressure measurement depends on matching the device cuff and size of your wrist. THE CUFF SHOULD NOT BE TOO SMALL OR TOO BIG.


(classification of the World Health Organization)

18 Repeated measurements are carried out at 5-minute intervals to recover the blood circulation and to avoid the injury due to blood flow interference. However, persons suffering from severe atherosclerosis, due to a significant loss of elasticity of blood vessels, need longer intervals between measurements (10-15 minutes).

This also concerns patients suffering from long-term diabetes. For more accurate determination of blood pressure it is recommended to carry out a series of three consecutive measurements and to calculate the average value of measurement results.

## POWER SUPPLY OF THE DEVICE

- 1 Remove the battery compartment. Do not use excessive force when removing the cover.
- 2 Install two "AAA". Make sure that polarity corresponds to signs (+) and (-) shown inside the compartment. Batteries are readily installed by pressing the end "-" on the spring.
- 3 Close the battery compartment cover. After the installation of batteries, set date and time.

 Battery Replacement Indicator

Replace all the batteries when the display continuously shows the battery replacement indicator " " or there is no indication on the display. The battery replacement indicator does not show a discharge degree.

Installation of new alkaline batteries provides up to 300 measurements. The supplied batteries are used for testing the device in the sale and service life may be shorter than that recommended alkaline batteries. Attention should be paid for the expiry date of the batteries especially in case the monitor is not used frequently.

Since neither the device nor the batteries are the waste that can be utilized at home, follow your national/local regulations for waste recycling and take them to corresponding collection facilities.

**ATTENTION!** Do not use rechargeable batteries or the batteries other than the ones with designated type. The specified performance of the product may not be retained if those batteries are used. Do not touch another person while touching the battery terminals of the monitor when changing batteries.

## SETTING DATE AND TIME


Date and time can be set after installing batteries. Setting the date and time guarantees the preservation of measurement results with indicated correct date and time.

If you do not want to enter date and time after batteries installation, press START button  five times.

### 1 Setting the Year

After batteries installation, the LCD display shows flashing digits of the year. Enter the current year using M1 button to increase the year value, and M2 button to decrease it. Press the START button  to confirm and to pass to the next step.

### 2 Setting the Month



Upon entering the current year, enter the current month (month digits flash on LCD display). Enter the current month using M1 button to increase the month value, and M2 button to decrease it. Press the START button  to confirm and to pass to the next step.

### 3 Setting the Day

Upon entering the current month, enter the current date (date digits flash on LCD display). Enter the current date using M1 button to increase the day value, and M2 button to decrease it. Press the START button  to confirm and to pass to the next step.


### 4 Setting the Time

Watch uses a 24-hour time format of day.

Upon entering the date, enter the hour using M1 and M2 buttons. Press the START button  to confirm the settings. Using M1 button to increase the minute value, and M2 button to decrease it. Press the START button  to confirm the settings. The clock has been set, and the unit is prepared for operation. LCD monitor will display the time even when the unit is switched off.

## CUFF PREPARATION

- 1 Holding your left hand with the palm facing up place a cuff on a wrist so that the meter display is on the palm side. Should you fail to put the cuff on your left wrist, place it on the right one.
- 2 Place the cuff on your hand so that its edge is spaced by 5-10 mm off the palm end. Place the meter on the center of your wrist.
- 3 Snugly wrap the cuff on your wrist without leaving any space between the cuff and the wrist. Make sure that the cuff comfortably fits your wrist.


Put the cuff on a bare wrist only. Be careful to prevent your cloth being caught by the cuff.

## CORRECT POSITION DURING MEASUREMENT

### Measurement in the sitting position

1. Sit down on a chair that you can be comfortably seated with your legs uncrossed, feet flat on the floor and back supported. Your arm should also be supported on the table.
2. Lightly raise your left hand with the palm facing up and put your elbow on the table.
3. Bring the cuff to the heart level placing your hand on a case or folded towel.

### In case of unavailability of a table

1. Sit down on a chair that you can be comfortably seated with your legs uncrossed, feet flat on the floor and back supported.
2. Bring the cuff to the heart level lightly pressing your left hand to the chest.
3. While measuring lightly support your left hand by the right one.

### Lie on your back

1. Lie on your back.
2. Bring the cuff to the heart level using a case and folded towel.

Try to use the same wrist and position for measurement.


If the cuff is positioned above/below the heart, any reading obtained will tend to be out of the allowance.

## MEASUREMENT PROCEDURE

Before measurement, take a few breaths and relax. During the measurement, do not talk and do not move.

### 1 Press the START «».

The unit will automatically start pumping the air into the cuff. When the cuff pressure reaches ~190 mm Hg, the unit will stop pumping the air.

Press the START «» button to stop forcibly the measurement: the device will stop inflation and quickly release the air and switch off.

The pressure in the cuff (the value shown on the display) starts to decrease, and heart rate is displayed in the form of icon ♥ (heart).

When the measurement is completed, the air is automatically released from the cuff.

The display shows systolic, diastolic blood pressure and pulse rate.

If irregular pulse rhythm is detected during the measurement, then, upon completion of the test, the irregular pulse rhythm icon will appear (two hearts, one of them flashing).

### 2 Press the M1 button or M2, button and the result will remain in the selected memory unit.

The number of the selected memory cell is shown on the display. The result is stored in the memory determined at the time the device is turned off. Only the blood pressure and heart rate icons will be saved, not the irregular pulse rhythm icon. If the measurement ended up with an error, the result will not be saved.

### 3 Press the START «Ⓢ» button to turn off the device.

If you forget to turn off the device, it will do so automatically after 3 minutes.

#### AUTOMATIC PRESSURE BOOSTING


If during measurement the initial cuff pressure (190 mm Hg) is not enough or arm is moved, the Device will stop measurement and pressurize the cuff to the next, higher pressure level. Automatic cuff pressure boosting will continue until measurement is completed successfully. This is not a malfunction.

To change over to manual pressure pump control, press START «Ⓢ» button immediately after beginning of air pumping into the cuff, and release it when the pressure reaches the level that is 50 mm Hg higher than the expected maximum systolic pressure. Air pumping will be stopped after releasing the button, if the pressure exceeds the value of 180 mm Hg. The pressure may be increased up to 300 mm Hg.

#### INDICATION OF IRREGULAR PULSE RHYTHM

The flashing symbol "♥", appeared on the display reports on an irregular pulse rhythm. The irregular pulse rhythm indicator occurrence can be also caused by body movements in the course of the measurement. If this symbol appears from time to time, refer to your doctor.

Do not perform several measurements in a row. This will cause numbing the hand and can affect the measurement result. Give your hand a break for at least 5 minutes.

## MEMORY FUNCTOIN

The measured values are automatically saved for later viewing in the selected memory unit. These memory units can be used to save separately the measurement results of two different persons or to save separately the results of morning and evening measurements.

Each unit can keep up to 30 measurement results and the average value thereof. When the number of kept values exceeds 30, the oldest readings are deleted to record new values.

Saved values are kept with the measurement date and time, if set. If you want to save time and date together with the measurement values, date and time should be set before the measurement.

When an error occurs (ERR), the measurement results are not saved.

## Viewing the Saved Data

1 Press M1 or M2 button to view the data saved.

2 Initially, the average value of results saved in the selected memory is displayed on the LCD screen. The «A» index stands for the average value.



If there is only one measurement result available, only this result is displayed, and if there are no results saved, the display shows the time.

The indication in the upper part of the display alternately changes from the number of memory data to the date and the time. Memory data will be shown for about 30 sec, and then the device automatically switches to the time display mode.

3 To view the next memory cell, press M1 button again (or M2, respectively, to view the content of memory M2). The result displayed in a cell No. 1 is the latest among the data saved in the selected memory.

4 Every press of M1 or M2 button results in a transition to the next memory cell. The bigger is the number of the memory cell, the older is the result.

5 Press the START «O» button to turn off the device.

## Deletion of Saved Data



The unit features the option of cleaning both a certain memory cell and the entire memory M1 (or M2) at once.

1 Select a value from the memory unit to be deleted or the average value (unit with index «A») to clear the whole memory unit.

2 Press and hold the button until the sign «←» appears on the display.

## INFORMATION ABOUT ERRORS

INDICATION	LIKELY CAUSE	METHODS OF CORRECTION
<p>Excess of pumping pressure</p> 	<p>Maximum allowable pressure: pressure can not be measured due to movement or talking during the measurement, in spite of the fact that the cuff was pumped at maximum.</p>	<p>Be calm and quiet during the measurement.</p>

<p>Measurement error</p> 	<p>During the measurement, a person was talking or moving.</p>	<p>Do not talk and do not move during the measurement.</p>
<p>Pumping error</p> 	<p>The cuff is put on incorrectly. The cuff has been damaged</p>	<p>Put on the cuff properly and repeat the measurement. If the Err-2 error repeats, that means the cuff has been damaged. Contact your local service center.</p>
<p>Air release error</p> 	<p>During the measurement, a person was talking or moving.</p>	<p>Turn off the device and perform the measurement again. Do not talk and do not move during the measurement.</p>
<p>Test mode</p> 	<p>You touched the START «» button when installing the batteries.</p>	<p>Turn off the device by pressing the START «» button and perform the measurement again.</p>
<p>No indication.</p>	<p>Batteries are discharged. Batteries are installed incorrectly. Connecting terminals are contaminated.</p>	<p>Replace all batteries with new ones. Install batteries properly. Wipe connecting terminals with a dry cloth.</p>
<p>Blood pressure is too low or too high.</p>	<p>The cuff is not at the heart level. The cuff has not been wrapped tightly around the wrist. Your shoulders and arms are strained. Probably your talking and moving your arm was the case at measurement.</p>	<p>Put on the cuff at the heart level. Attach the cuff correctly.  Relax and repeat the procedure.  Stay calm and be silent during measurement.</p>
<p>Blood pressure values vary between measurements. Measurement values are too low (high).</p>	<p>Blood pressure values continuously change depending on the measurement time and nervous excitation.</p>	<p>Take some deep breaths to relax before the measurement.</p>

If, despite the above-given recommendations, you fail to obtain the right measurements, stop the operation and contact the service center (addresses and telephone numbers of authorized organizations are provided in the warranty certificate). Do not attempt to adjust the device internal mechanism on your own.

## TECHNICAL SPECIFICATIONS

Operating Principle	Oscillometric method
Indicator	4-line LCD display
Indicating Range: cuff pressure, mmHg	0-300
Measuring Range: cuff pressure, mmHg pulse rate, bpm	50-250 (systolic), 40-180 (diastolic) 40 - 160
Accuracy: cuff pressure, mmHg pulse rate, %	$\pm 3$ $\pm 5$
Inflation	Automatic inflation (air pump)
Deflation	Automatic Deflation (electric valve)
Electric supply voltage, V	3
Electric supply type	Two AAA (LR03) batteries
Power Consumption, W	2
Memory	2 units, 30 values per each + average value
Operating Condition temperature, °C humidity, % Rh	+10 to +40 85 and lower
Storage Condition temperature, °C humidity, % Rh	-5 to +50 85 and lower
Cuff	Cuff WS-820
Cuff Size	adult (circumference of the wrist 12,5-21,5 cm)
Overall dimensions: Size (electronic unit), mm Weight (without packaging, bags and batteries), g	63,5 x 71 x 33 100
Manufacture year:	specified on the device housing in the device serial number after the letters «SN»
Protection class IP:	IP20; Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.
Protection against electric shock:	Internally powered equipment, Type BF applied part (cuff)
Mode of operation:	Continuous operation
Classification:	Internally Powered Equipment

### KEY TO SYMBOLS:

Important: Please read the Manual

European Authorized Representative

Type BF applied part

Name and address of the manufacturer

Environmentally safe packaging

When utilizing the waste, refer to current rules applicable in your region

0123 Compliance with the Directive 93/42/EEC

Keep dry

The measurement accuracy of the product has been proven according to ISO 81060-2 protocol. In the clinical study, K5 was used for the determination of diastolic pressure values at all auscultatory measurements.

This product is intended to be used in the condition of one atmospheric pressure.

Revision date of the present Manual is indicated on the last page as EXXX/YYMM/XX, where YY is the year and MM is the month of revision.

## WARRANTY

1. The manufacturer guarantees the warranty period of 5 years for the device from the sale date provided that the consumer observes operation, transportation and storage requirements. The warranty period for the cuff and the power source is 12 months from the sale date.
2. Warranty liabilities are documented with the warranty certificate upon selling the device to the buyer. The guarantee is valid provided that the device has not been opened or damaged by the buyer.
3. Addresses of organizations engaged in the warranty service are specified in the warranty certificate.

## CARE, STORAGE, REPAIR AND DISPOSAL

- 1 This device should be protected from excessive temperature or moisture, direct sunlight, strokes, dust, lint and vibration. **THE DEVICE IS NOT WATERPROOF!**
- 2 Do not keep or do not use the device in close proximity to heaters and open flame.
- 3 If the device has not been used for a long time, remove the batteries. Leaking of batteries can cause damage to the device and terminate the warranty. **KEEP BATTERIES AWAY FROM CHILDREN!**
- 4 Do not contaminate the device and protect it from dust. Check and ensure the monitor is not dirty before use. In case the cleaning is necessary, clean the body with a cloth moistened with water and a mild detergent, and then wipe dry.
- 5 Do not allow the contact between the device and its parts with water, solvents, alcohol, and gasoline.
- 6 Keep the cuff away from sharp objects, and do not try to pull out the cuff.
- 7 Do not expose the device to strong strokes and do not throw it.
- 8 The device does not contain any adjustment controls for settings. Unauthorized opening of the electronic device is forbidden. If needed, repair the device only in specialized organizations.
- 9 On the expiry of the specified operation term, refer to specialists (specialized repair organizations) on a periodic basis to check the technical condition of the device. We suggest that you have your monitor checked every 2 years.
- 10 When utilizing the waste, refer to current rules applicable in your region. No special utilization conditions are specified by the manufacturer for this device.
- 11 The cuff is resistant to repeated sanitation. The cuff internal fabric surface (being in contact with a patient's wrist) can be treated with a cotton swab moistened in a 3% solution of hydrogen peroxide. Partial discoloration of the cuff covering tissue is possible if used for a long time. Do not wash the cuff and do not treat it with a hot iron.
- 12 Make sure that the monitor is stored in a secure place where children or pets cannot reach.
- 13 Do not replace the parts, disassemble or modify the monitor as it causes malfunction of the monitor or danger to users.
- 14 Stop using the monitor immediately and contact your dealer or the manufacturer in case any visible damage is found on the monitor.

## CERTIFICATION AND STATE REGISTRATION

The production of devices is certified pursuant to international standards such as ISO 9001, ISO 13485, ISO 14001. The device meets international standards IEC60601-1+A1+A2, IEC60601-1-2+A1/IEC60601-1-2 CISPR 11+A1. Device comply with the requirements of European Directive MDD 93/42/EEC.

Produced by Nihon Seimitsu Sokki Co., Ltd.

Address: 2508-13 Nakago Shibukawa Gunma 377-0293 Japan

website: <http://www.nissei-kk.co.jp/english/>

EC-Representative: Nissei Healthcare (UK) Ltd.

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website: <http://www.nisseihealthcare.com>

## TECHNICAL DESCRIPTION

WS-820 complies with the EMC, electromagnetic compatibility, standard, IEC60601-1-2. Refer to the following tables for specific information regarding compliance to the standard. WS-820, as a medical electrical equipment, needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

- Portable and mobile RF communications equipments can affect the device.
- The use of accessories other than those specified in this manual may result in increased emissions or decreased immunity of the device.
- WS-820 should not be used adjacent to or stacked with other equipment.

Table 1 - Guidance and manufacturer's declaration - electromagnetic emissions -

WS-820 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-820 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	WS-820 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  Note: Operate by DC Power
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Table 2 - Guidance and manufacturer's declaration - electromagnetic immunity -

WS-820 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-820 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	Note: Operate by DC Power.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	Note: Operate by DC Power.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	N/A	Note: Operate by DC Power.
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.			

Table 4 - Guidance and manufacturer's declaration - electromagnetic immunity -

WS-820 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-820 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Portable and mobile RF communications equipment should be used no closer to any part of WS-820, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms, 150 kHz to 80 MHz 3 V/m, 80 MHz to 2.5 GHz	N/A  3 V/m	Recommended separation distance $d=[3.5/V1]\sqrt{P}$  $d=[3.5/E1]\sqrt{P}$ , 80MHz-800MHz $d=[7/E1]\sqrt{P}$ , 800MHz-2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<sup>a</sup> Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which WS-820 is used exceeds the applicable RF compliance level above, WS-820 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating WS-820.			
<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 6 - Recommended separation distances between portable and mobile RF communications equipment and WS-820 -

WS-820 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of WS-820 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and WS-820 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz to 80 MHz, $d=[3.5/V1]\sqrt{P}$	80 MHz to 800 MHz, $d=[3.5/E1]\sqrt{P}$	800 MHz to 2.5 GHz, $d=[7/E1]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

CE 0123

**NIHON SEIMITSU SOKKI CO., LTD.**

2508-13 Nakago Shibukawa Gunma 377-0293 Japan

Nissei Healthcare (UK) Ltd.

Rede House, New Barn Lane Henfield West Sussex BN5 9SJ UNITED KINGDOM

® Registered trademark.

© Copyright 2017.

E253/1701/02

A111623-1\_6

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Produsul companiei Nihon Seimitsu Sokki Co., Ltd.